

VAV SYSTEMS (SAFE AND SECURE)

PLASMA EXTRACTOR



Plasma Extractor is designed to Extract Blood component from centrifuge Bags. The spring loaded front panel applies pressure on the collected bag causing the liquid to contain in transfer bag. Its sturdy construction makes it durable, easy to use and portable, the unit is available in both version automatic and manual.

Automatic Plasma Expressor. The spring loaded front panel applies pressure on the collected bag causing the liquid to contain in transfer bag. These have an attractive front panel bezel and caters to the requirements of hospitals, research clinics etc.

Features:

- IR Sensor
- Motor activated clamping
- Audiovisual Alarm
- Spring Loaded Acrylic Plate
- Uniform Pressure
- Automatic Control
- POWER 230V AC

SR.NO.	SIZE	WEIGHT	ACRALIC THICKNESS./SIZE
1	6/10/8 inch(L/W/H) (MANUAL)(VSPX50)	4 KG.	10MM./6/8 INCH
2	6/10/8 inch(L/W/H) (AUTOMETIC)(VSPX 100)	4.5KG.	10MM./6/8 INCH