

VAV INSTRUMENTS (SAFE AND SECURE)



QUICK FREEZERS

Our high-performance -35°C freezers provide long-term storage for specialty plasma products, delivering protection against interior warm-up that can result from heavy freezer use. Our -40°C plasma freezers maintain safe, uniform storage conditions even at ambient temperatures of up to +45°C.

- -35°C and -40°C units with storage capacities from 80 liter to 600 liter.
- Microprocessor based Digital display Independent bottom temperature control with audio visual alarm for high and low temperature. And Power failure alarm.
- CFC free refrigeration system.
- Start up delay protection.
- High intensity insulation foam.
- Lockable casters- easy to move.
- It has washable filter.
 - Positive latch with key lock security; prevents door 'pop-open' on two-door units.
 - Anti-sweat door frame heaters with magnetic compression gaskets to provide a tight door seal.
- dual door for zero cooling loss.

MODELS	VIPF80	VIFP140	VIPF250	VIPF300	VIPF450
TEMPERATURE	-35°C	-35°C	-35°C	-35°C	-35°C
NOMBER OF TRAYS	2	4	5	6	7
SIZE(W*D*H)cms Inner size.(app.)	44*40*61	49*50*82	57*49*119	60*55*120	63*65*142
VOLUME	150ltr.	250 ltr.	300 ltr.	450 ltr.	600 ltr.

SIZE MAY CHANGE WITHOUT PRIOR NOTICE.